University Center A. Elhafidh Boussouf Mila Institute of Science and Technology 1st year ST 2023/2024 Chemistry Experiments

Experiment 2: Preparation of Solutions in Chemistry

### Introduction

**Solution :** A solution is a homogeneous mixture of two or more pure substances that do not interact with each other. A solution is formed when a substance, called a solute, is dissolved in another substance, called the solvent.

Solute : Substance which dissolves in a solution.

Solvent : Substance which dissolves another to form a solution (water).

Solution = Solute + Solvent

#### Ways of expressing concentration

There are numerous ways to describe the concentration of a solution, and they are a useful ways to describe solutions concentrations in chemistry reactions.

Molarity	$C_{\rm M} = \frac{n_{solute}}{V_{solution}} \; (\frac{mol}{L})$
Mass concentration	$T = \frac{m_{solute}}{V_{solution}} \left(\frac{g}{L}\right)$
Normality	$C_{\rm N} = \frac{n_{eq.g}}{V_{solution}} \left(\frac{eq.g}{L}\right)$
Molality	$C_{\rm m} = \frac{n_{solute}}{m_{solvent}} \left(\frac{mol}{Kg}\right)$

### **Relation between Normality and Molarity**

Normality and Molarity are two important and commonly used expressions in chemistry. They are used to indicate the quantitative measurement of a substance. But what relation does Molarity have with Normality ?

$$C_N = Z C_M$$

Where Z is a constant with a number of states

- For acid : Z is the number of H<sup>+</sup> protons that the acid can lose. (HCI : Z = 1) (H<sub>2</sub>SO<sub>4</sub>: Z = 2).
- For base : Z is the number of hydroxides OH<sup>-</sup> that the base can lose (NaOH : Z=1 ; BaOH<sub>2</sub> : Z=2)
- **For oxidation and reduction :** Z is the number of transferred electrons.  $[MnO_4^- + 5e_7 + 8OH_7 Mn^{2+} + 4H_2O] : Z = 5$
- **For salts :** Z is the number of metal atoms in its valence.  $[Al_2(SO_4)_3 2Al_3 + + 3SO_4^{2-}] : Z = 2 \times (+3) = 6$

#### **Objective of the experiment**

- 1. Recognising and using the equipment and tools for preparing solutions.
- 2. How to prepare a solution from sodium hydroxide (NaOH) by dissolving.
- 3. How to dilute Hydrochloric acid (HCI) solution.

#### **Materials and Chemicals**

Materials	Chemicals
Graduated pipette or graduated cylinder	Distilled water
Volumetric flask	• Hydrochloric acid (HCI)
• Spatula - Watch Glass - Funnel	• Sodium hydroxide (NaOH)
Analytical balance	

#### Procedure

**1. Preparation of solution from solid :** This is how to make a chemical solution using a solid dissolved in a liquid. We need to calculate the mass of the solid required using with the equation :  $\mathbf{m} = \mathbf{C}_{\mathbf{M}} \cdot \mathbf{M} \cdot \mathbf{V}$ 

**Question :** Calculate and describe how to prepare 100 ml of sodium hydroxide (NaOH) solution with a molar concentration of 0.1 mol/L ( $M_{NaOH} = 40 \text{ g/mol}$ )?


- Calculate the mass of sodium hydroxide (NaOH) needed to prepare 100 ml of 0.1 mol/L NaOH solution.

$$C_{M} = \frac{n}{V} = \frac{m}{M.V} \rightarrow m = C_{M}.M.V$$

$$m_{NaOH}(g) = C_{\text{Integer}}(\frac{mol}{L}) \cdot M_{NaOH}(\frac{g}{mol}) \cdot V_{\text{Integer}}(L)$$
$$m_{NaOH}(g) = 0.1 \times 40 \times 0.1 = 0.4 \text{ g}$$

- Weight the mass of 0.4 g of NaOH.
- Fill a clean 100 mL volumetric flask at third of it with water.

- Add 0.4g of NaOH to this volumetric flask using a funnel.
- Stir the mixture until the NaOH is completely dissolved.
- Completely fill the volumetric flask with distilled water to the measuring line.
- Close the volumetric flask, then mix to obtain homogeneous solution of NaOH.



Scheme 1. How to prepare solutions from solid

**2. Preparation of solution from liquid :** This is how to make a chemical solution through dilution method. We need to calculate the volume of the starting solution required using with the equation :

$$C_i . V_i = C_f . V_f \rightarrow V_i = \frac{C_f . V_f}{C_i}$$

**Question :** Calculate and describe how to prepare 100 ml of 0.1 mol/L hydrochloric acid HCI from concentarted HCl solution ?

- Read the information on the concentrated HCl(commercial); the density is 1.18, the purity rate is 37%, and the molar mass is 36.5 g/mol.
- Calculate the mass of the concentrated HCl.

We have d= 1.18 so  $\rho=1.18$  Kg/L i.e. one L of HCI weighs 1180 g.

- Determine the amount of pure concentrated HCl in grammes.

$$m_{HCl}$$
(Concentrated) =  $\frac{1180 \times 37}{100}$  = 436.6 (g)

- Determining the concentration of concentrated HCl.


$$C_{HCl} = \frac{n_{HCl}}{V} = \frac{m_{HCl}}{M_{HCl} \cdot V} = \frac{436.6}{36.6 \times 1} = 11.96 \left(\frac{\text{mol}}{\text{L}}\right)$$

- Calculating the volume required of concentrated HCl to prepre the needed solution

$$V_{HCl} = \frac{C_{\rm f} \cdot V_{\rm f}}{C_{\rm HCl}} = \frac{0.1 \times 100}{11.96} = 0.83 \,({\rm ml})$$

- Fill a clean 100 mL volumetric flask at third of it with water.

- Take out 0.83 mL of concentrated HCl using the graduated pipette.
- Transfer it to the volumetric flask.
- Completely fill the volumetric flask with distilled water to the measuring line.
- Close the volumetric flask, then mix to obtain homogeneous solution of HCl.



Scheme 2. How to prepare solutions through dilution

## Answer the questions

Calculate and describe how to prepare the following solutions :

- 1. 100 mL of 0.1 M KOH using solid KOH.
- 2. 200 ml of 0.25N (eq.g/L)  $H_2SO_4$  using concentrated  $H_2SO_4$
- 3. 250 ml of a 0.3 M CH<sub>3</sub>COOH from 0.5 M CH<sub>3</sub>COOH.

# We give

KOH: (M = 56.11 g/mol - % = 85)H<sub>2</sub>SO<sub>4</sub>: (M = 98 g/mol - d = 1,18 - % = 96)CH<sub>3</sub>COOH: (M = 60 g/mol - d = 1.05 - % = 99)