

Practical work N°2. Preparation of aqueous solutions in chemistry

Introduction

Solution: A solution is a homogeneous mixture of two or more pure substances that do not interact with each other. A solution is formed when a substance, called a solute, is dissolved in another substance, called the solvent.

Solute: Substance which dissolves in a solution.

Solvent: Substance which dissolves another to form a solution (water)

$$\text{Solution} = \text{Solute} + \text{Solvent}$$

Ways of expressing concentration

There are numerous ways to describe the concentration of a solution, and they are a useful ways to describe solutions concentrations in chemistry reactions.

Molarity	$C_M = \frac{n_{\text{solute}}}{V_{\text{solution}}} \left(\frac{\text{mol}}{L} \right)$
Mass concentration	$T = \frac{m_{\text{solute}}}{V_{\text{solution}}} \left(\frac{g}{L} \right)$
Normality	$C_N = \frac{n_{\text{eq.g}}}{V_{\text{solution}}} \left(\frac{\text{eq.g}}{L} \right)$
Molality	$C_m = \frac{n_{\text{solute}}}{m_{\text{solvent}}} \left(\frac{\text{mol}}{Kg} \right)$

Relation between Normality and Molarity

Normality and Molarity are two important and commonly used expressions in chemistry. Both are used to indicate the quantitative measurement of a substance. But what relation does Molarity have with Normality?

$$C_N = Z C_M$$

Where Z is a constant with a number of states

✚ **For acid:** Z is the number of H^+ protons that the acid can lose. (HCl : Z = 1) - (H_2SO_4 : Z = 2).

✚ **For base :** Z is the number of hydroxides – OH^- that the base can lose (NaOH : Z=1 ; $Ba(OH)_2$: Z=2)

✚ **For oxidation and reduction:** Z is the number of transferred electrons.
[$MnO_4^- + 5e^- + 8OH^- \rightarrow Mn^{2+} + 4H_2O$]: Z = 5

✚ **For salts:** Z is the number of metal atoms in its valence.
[$Al_2(SO_4)_3 \rightarrow 2Al^{3+} + 3SO_4^{2-}$] : Z = $2 \times (+3) = 6$

Objective of the experiment

1. To identify and use the equipment and tools required for preparing solutions.
2. To prepare a sodium hydroxide (NaOH) solution by dissolution.
3. To dilute a hydrochloric acid (HCl) solution.

Materials and Chemicals

Materials	Chemicals
• Graduated pipette or graduated cylinder • Volumetric flask • Spatula - Watch Glass - Funnel • Analytical balance	• Distilled water • Hydrochloric acid (HCl) • Sodium hydroxide (NaOH)

Procedure

1. Preparation of solution from a solid:

This section describes how to prepare a chemical solution using a solid dissolved in a liquid.

We need to calculate the mass of the solid required using the equation: $m = C_M \cdot M \cdot V$

Question: Calculate and describe how to prepare 100 ml of sodium hydroxide (NaOH) solution with a molar concentration of 0.1 mol/L ($M_{NaOH} = 40 \text{ g/mol}$)?

- Calculate the mass of sodium hydroxide (NaOH) needed to prepare 100 ml of 0.1 mol/L NaOH solution.


$$C_M = \frac{n}{V} = \frac{m}{M \cdot V} \rightarrow m = C_M \cdot M \cdot V$$

$$m_{NaOH} (g) = C_{المحلول} \left(\frac{mol}{L} \right) \cdot M_{NaOH} \left(\frac{g}{mol} \right) \cdot V_{المحلول} (L)$$

$$m_{NaOH} (g) = 0,1 \times 40 \times 0.1 = 0.4 \text{ g}$$

- Weight the mass of 0.4 g of NaOH.
- Fill a clean 100 mL volumetric flask at third of it with water.

- Add 0.4g of NaOH to this volumetric flask using a funnel.
- Stir the mixture until the NaOH is completely dissolved.
- Completely fill the volumetric flask with distilled water to the measuring line.
- Close the volumetric flask, then mix to obtain homogeneous solution of NaOH.



Scheme 1. How to prepare solutions from solid

2. Preparation of solution from a liquid:

This section explains how to prepare a chemical solution through dilution method. We need to calculate the volume of the starting solution required using with the equation :

$$C_i \cdot V_i = C_f \cdot V_f \rightarrow V_i = \frac{C_f \cdot V_f}{C_i}$$

Question: Calculate and describe how to prepare 100 ml of 0.1 mol/L hydrochloric acid HCl from concentrated HCl solution ?

- Read the information on the concentrated HCl (commercial); the density is 1.18, the purity rate is 37%, and the molar mass is 36.5 g/mol.
- Calculate the mass of the concentrated HCl.

We have $d = 1.18$ so $\rho = 1.18 \text{ Kg/L}$ **i.e.** one L of HCl weighs 1180 g.

- Determine the amount of pure concentrated HCl in grammes.

$$m_{\text{HCl}}(\text{Concentrated}) = \frac{1180 \times 37}{100} = 436.6 \text{ (g)}$$

- Determining the concentration of concentrated HCl.


$$C_{\text{HCl}} = \frac{n_{\text{HCl}}}{V} = \frac{m_{\text{HCl}}}{M_{\text{HCl}} \cdot V} = \frac{436.6}{36.6 \times 1} = 11.96 \left(\frac{\text{mol}}{\text{L}} \right)$$

- Calculating the volume required of concentrated HCl to prepare the needed solution

$$V_{\text{HCl}} = \frac{C_f \cdot V_f}{C_{\text{HCl}}} = \frac{0.1 \times 100}{11.96} = 0.83 \text{ (ml)}$$

- Fill a clean 100 mL volumetric flask at third of it with water.

- Take out 0.83 mL of concentrated HCl using the graduated pipette.
- Transfer it to the volumetric flask.
- Completely fill the volumetric flask with distilled water to the measuring line.
- Close the volumetric flask, then mix to obtain homogeneous solution of HCl.



Scheme 2. How to prepare solutions through dilution

Given

KOH: ($M = 56.11 \text{ g/mol}$ - 85%)

H_2SO_4 : ($M = 98 \text{ g/mol}$ - $d = 1.18$ - 96%)

CH_3COOH : ($M = 60 \text{ g/mol}$ - $d = 1.05$)